


pLenti-Puro+NAMPT • Flag plasmid

Cat No.: VGP-10135-001P2

Gene/insert name	Homo sapiens nicotinamide phosphoribosyltransferase
Official Symbol	NAMPT
Alt name	1110035O14Rik; PBEF; PBEF1; VF; VISFATIN
Insert size	1476 bp
Gene ID	10135
Accessions	NM_005746.3
Tag	Flag
Backbone size	8038 bp
Cloning site 5'	Xba I
Cloning site 3'	Not I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

