


pLenti -GFP-Puro+RAMP1 plasmid

Cat No.: VGP-10267-71H4

Gene/insert name	Homo sapiens receptor activity modifying protein 1, transcript variant 1, mRNA
Official Symbol	RAMP1
Alt name	N/A
Insert size	922 bp
Gene ID	10267
Accessions	NM_005855.4
Tag	N/A
Backbone size	8177 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

