


pLenti-GFP-Puro+LINC01843 plasmid

Cat No.: VGP-101927934-001P1

Gene/insert name	Homo sapiens long intergenic non-protein coding RNA 1843, long non-coding RNA
Official Symbol	LINC01843
Alt name	Lnc-APUE
Insert size	1108 bp
Gene ID	101927934
Accessions	NR_105045.1
Tag	N/A
Backbone size	7981 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

