


pLenti-GFP-Puro+ELFN1-AS1 plasmid

Cat No.: VGP-101927125-001P1

Gene/insert name	Homo sapiens ELFN1 antisense RNA 1, transcript variant 1, long non-coding RNA
Official Symbol	ELFN1-AS1
Alt name	MYCLO-2
Insert size	1008 bp
Gene ID	101927125
Accessions	NR_120508.1
Tag	N/A
Backbone size	8081 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	GFP

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

