


pLenti -GFP-Puro+ABCC4 • Flag plasmid

Cat No.: VGP-10257-001P1

| | |
|--------------------------------|--|
| Gene/insert name | Homo sapiens ATP-binding cassette, sub-family C , member 4, mRNA |
| Official Symbol | ABCC4 |
| Alt name | MOAT-B; MOAT-B; MRP4 |
| Insert size | 2580 bp |
| Gene ID | 10257 |
| Accessions | BC041560.1 |
| Tag | Flag |
| Backbone size | 8171 bp |
| Cloning site 5' | Xba I |
| Cloning site 3' | BamH I |
| Growth strain(s): | Stbl 3 |
| High or low copy | High |
| Bacterial resistance(s) | Ampicillin |
| Growth temperature | 37 °C |
| Vector type | Lentiviral, Mammalian Expression |
| Selectable markers | Puromycin |
| Fluorescent protein | GFP |

Map



| | |
|---------------------------------|---|
| Buffer | TE buffer |
| Storage | Store at -20°C up to 24 months |
| The plasmid is ready for | Transfection |
| Note | For research use only, not for use in diagnostic procedure. |

