


pLenti-Puro+SP140(h) • 3xFlag plasmid

Cat No.: VGP-11262-36H8

Gene/insert name	Homo sapiens SP140 nuclear body protein,transcript variant 1
Official Symbol	SP140
Alt name	LYSP100; LYSP100-A; LYSP100-B
Insert size	2604 bp
Gene ID	11262
Accessions	NM_007237.5
Tag	3xFlag
Backbone size	7428 bp
Cloning site 5'	Xba I
Cloning site 3'	BamH I
Growth strain(s):	Stbl 3
High or low copy	High
Bacterial resistance(s)	Ampicillin
Growth temperature	37 °C
Vector type	Lentiviral, Mammalian Expression
Selectable markers	Puromycin
Fluorescent protein	W/O

Map



Buffer	TE buffer
Storage	Store at -20°C up to 24 months
The plasmid is ready for	Transfection
Note	For research use only, not for use in diagnostic procedure.

